1. Warnings and cautions

WARNINGS

- Insufficient, reduced or lack of ventilation may lead to brain damage in the ventilated patient.
- Should only be used by personnel trained in performing resuscitation. It is important to obtain a tight seal between the face and the mask. Otherwise, this may cause reduced ventilation or no ventilation. Make sure that such personnel are made familiar with the contents of these instructions for use.
- Watch the movement of the chest, and listen for the expiratory flow from the valve to check the efficiency of the ventilation. Failure to do so may result in insufficient ventilation.
- Switch IMMEDIATELY to mouth-to-nose or mouth-to-mouth ventilation in accordance with local guidelines if efficient ventilation cannot be obtained with the device. Failure to do so may result in insufficient ventilation.
- Do not smoke or use open flames when using oxygen - fire may result.
- Do not use the resuscitator in toxic or hazardous atmospheres as this may harm the patient. See Section 5.4 for further details.
- Once the override cap has been mounted on the pressure-limiting valve and the whole hand is used for compressing the Ambu Mark IV Baby resuscitator, the pressure-limiting outer cover is prevented from expanding, meaning that in case of extremely heavy compression, a pressure as high as 8-10 kPa (80-100 cmH₂O) is created in patients with low lung compliance.

CAUTIONS

- According to US federal regulations, this product may only be sold by or on the order of a physician (only applies in the USA and Canada).
- It is important to check that all parts of the resuscitator are in good condition, that it is properly assembled, and that a complete functional test is carried out every time the device has been disassembled - also before the initial use.
- This functional test, which is described in chapter 8.1, should also be carried out every time the device is taken into use. Check of correct function of the pressure-limiting valve according to Section 8.2 should also be performed immediately before use.
- If the resuscitator and its accessories are kept on standby for emergency situations, the entire equipment should be checked regularly to ensure that it is functional.
- Avoid the use of or contact with oil or grease as they may damage the materials of the resuscitator. In case of contact, the surface must be thoroughly cleaned.
- Do not use oil or grease in close proximity to oxygen equipment - fire may result.
- Do not force the outer cover of the bag to bulge locally by closing the patient connector and pressing the bag so hard with your hand that a local bulge appears on the outer cover - this has an extreme impact on the material.
- The tube connector may only be used for pressure monitoring. The cap must always be put on the connector when pressure monitoring does not take place.

2. Application

The **Ambu Mark IV Baby resuscitator**, is intended for pulmonary resuscitation of neonates, infants and children with a body weight of up to 20 kg, approx. 4-5 years.

3. Specifications

The **Ambu Mark IV Baby resuscitator** complies with the following standards:

ASTM F 920-93 and EN ISO 10651-4:2002.

The Ambu Mark IV Baby resuscitator is in conformity with Directive 93/42/EEC concerning medical devices.

Application: Neonates, infants and children with a body weight of up to 20 kg, approx. 4-5 years.

Max. volume delivered: Approx. 300 ml

Length x bag diameter with patient valve: 265 mm x 85 mm

Weight of resuscitator with patient valve (without mask):
190 g

Max. frequency of ventilation: Depends on the tidal volume used. The expansion rate of the bag is more than sufficient for all practically applicable frequencies.

Patient connector: 22/15 mm (ISO)

Expiratory connector (for PEEP valve connection):
30 mm male (Ambu specified)

Dead volume: < 6 ml

Leaks (both ways): Non-measurable

Inspiratory resistance (with or without oxygen reservoir tube):

At 50 l/min: approx. -0.37 kPa (-3.7 cmH₂O)

At 5 l/min: approx. -0.05 kPa (-0.5 cmH₂O)

Pressure-limiting valve:

Opening pressure: approx. 4 kPa (40 cmH₂O)

Maximum attainable delivery pressure with compliance of 0.01 l/kPa when the pressure limiting valve is overridden:
approx. 11.5 kPa (115 cmH₂O).

Expiratory resistance:

At 50 l/min: approx. 0.45 kPa (4.5 cmH₂O)

At 5 l/min: approx. 0.06 kPa (0.6 cmH₂O)

Oxygen reservoir tube:

Length: 250 mm

Volume: 100 ml

Oxygen reservoir bag:

Inspiratory resistance:

With 50 l/min: approx. -0.1 kPa (-1.0 cmH₂O)

With 5 l/min.: approx. -0.01 kPa (-0.1 cmH₂O)

Volume: approx. 1500 ml

Application temperature: -18°C to 50°C (0°F to 122°F)

Storage temperature: -40°C to 70°C (-40°F to 158°F)

Parts / Materials

Inner bag	Silicone rubber
Outer cover	Silicone rubber
Patient valve housing	PSU
O-ring (internal part of patient connector)	Silicone rubber
Patientconnector	PSU
Pressure-limiting valve	PSU/stainless steel
Override cap	Silicone rubber
Expirations connector	Reinforced PP
Cap for manometer port	PP/EPDM rubber

Inlet valve housing	POM
Screw cap (two parts)	POM
Reservoir tube	PE
Valve discs	Silicone rubber
Ambu mask shell	PSU
Ambu mask 0A edge	Chloroprene with sealed foam inner lining
Ambu mask 0 edge	Silicone rubber
Rendell-Baker mask	Rubber
Hanging strap	Silicone rubber
Nameplate	Aluminium
Window with snap fastener	PSU
Oxygen reservoir bag (optional accessory)	PES/TPU
Reservoir housing incl. cap and screw cap (optional accessory)	POM/Reinforced PP
Adapter	Silicone rubber
Connector	POM

4. Operating principles

① shows how the ventilation air flows into the resuscitator and to and from the patient during manual ventilation as well as when the patient respire spontaneously. ① shows the operating principle of using an open reservoir tube 1.7 and oxygen reservoir bag (optional) 1.8.

O₂ reservoir system: During compression of the resuscitator, oxygen flows through the O₂ connector into the reservoir. When the resuscitator is released, oxygen is drawn from the reservoir into the resuscitator itself.

① 1.1: Pressure-limiting valve, 1.2: Air, 1.3: Excess oxygen, 1.4: Oxygen in, 1.5: Patient, 1.6: Expiratory.

The valve disc guide pin in the patient valve is a visible indicator of correct functioning of the unit. The user should always - before as well as during use - check that the valve guide moves forward on inspiration and backward on expiration. See ②.

② 2.1: Guide pin (forward), 2.2: Insufflation or spontaneous inspiration, 2.3: Guide pin (backward), 2.4: Expiration

5. Instructions for use

5.1. Ventilation with face mask

Clean mouth and airways using recommended techniques.

Use approved techniques for placing the patient's head, obtaining the correct head position for open airways. Place the mask so that it makes an airtight fit with the patient's face. Hold the mask firmly against the face while maintaining correct tilting of the head backward for the airways to be kept open.

WARNING

It is very important to practice correct application of the face mask before any attempt is made to use the resuscitator. Otherwise, the result may be no or reduced airflow.

Squeeze the bag with the tip of your fingers.

While insufflating, check that the patient's chest rises. See ③.

Loosen the grip on the bag. Listen for the expiratory air from the patient valve, and watch the chest lowering. See ③. Check at the same time whether the valve disc guide pins are moving in the transparent valve housing. If you feel constant resistance against the insufflation, check for blocked airways, or check the head's backward-tilted position.

The correct ventilation frequency may vary. Follow the current nationally or internationally recommended ventilation frequencies.

A pressure gauge can be connected to the tube connector on the side of the patient valve for monitoring the ventilation pressure. Remove the cap and connect the tube from the pressure gauge. See ⑤.

WARNING 

The tube connector may only be used for pressure monitoring. The cap must always be put on the connector when pressure monitoring does not take place.

In case of the patient vomiting under mask ventilation, the patient's airways must be immediately cleared of vomit. Squeeze the bag freely a few times without the mask being in contact with the face. Check that insufflation and expiration are unhindered.

The patient valve can be disassembled and cleaned if large volumes of vomit block the free airflow in the valve.

Perform functional test, see Section 8.1.

Administer oxygen according to medical indications. See Sections 5.3 and 5.4.

5.2. Pressure limitation - see ⑦

The patient valve is fitted with a pressure-limiting valve, which opens when the pressure in the valve reaches approx. 4.0 kPa (40 cmH₂O). This happens if the neonate's lungs are excessively full, the insufflation is too powerful or the airways are blocked. When the pressure-limiting valve opens, excess ventilation air is released into the atmosphere.

If a pressure above 4.0 kPa (40 cmH₂O) is required, the pressure-limiting valve can be overridden by pressing the override cap 7.3

onto the pressure-limiting valve. 7.1 shows the patient valve with the override cap off. 7.2 shows the patient valve with the override cap on. Mark IV Baby is supplied with the closing cap off.

The insufflation pressure can be fairly precisely controlled via the number of fingers used to compress the bag, see (4).

By using a tracheal catheter immersed in water you can get a sense of the feel of the insufflation pressure through the thin outer cover, see ⑥.

WARNING 

The pressure-limiting valve must never be overridden, unless a medical and professional evaluation indicates the necessity. High ventilation pressures may cause lung rupture in certain patients.

5.3. Oxygen enrichment - see ⑧

Administer oxygen according to medical indications.

- ⑧ 8.1: Oxygen connector (6 mm - this opening may only be occluded if a closed reservoir is mounted using an adapter),
- 8.2: Inlet connector (this opening may never be occluded),
- 8.3: Oxygen reservoir tube, 8.4: Adapter,
- 8.5: Oxygen reservoir bag, 8.6: Oxygen connector on oxygen reservoir bag (6 mm).

Oxygen enrichment is performed by connecting the oxygen supply tube to the connector on the back of the inlet valve and connecting the oxygen reservoir tube to the inlet connector.

When using the oxygen reservoir bag, use an adapter. Place the adapter over the inlet connector and close the oxygen connector using the cap (see ⑨). Oxygen supply takes place through the oxygen connector on the reservoir housing. Obtainable oxygen concentrations using an oxygen reservoir bag are as minimum what can be obtained by using the oxygen reservoir tube, see ⑪).

Oxygen enrichment of the inspiratory air in case of spontaneous respiration

- Without an oxygen reservoir tube:
Attach a mask to the resuscitator, set the oxygen flow to between 2 and 4 l/min, and place the mask over the child's mouth and nose, creating an airtight seal. If the spontaneous respiration is sufficient, an oxygen concentration close to 100% will be obtained.
- With an oxygen reservoir tube:
The end of the resuscitator with the oxygen reservoir tube that normally faces away from the patient can also be used for oxygen enrichment of the inspiratory air in case of spontaneous respiration. Again, the oxygen flow should be set to between 2 and 4 l/min, and the free end of the oxygen reservoir tube (1.3) should be held as tightly as possible against the child's head.

Again, the oxygen concentration will be very close to 100% at the end of the tube. Depending on the distance between the reservoir tube and the child's head, there will be significantly more or less inhalation of atmospheric air, and thus lower oxygen concentrations should be expected.

5.4. Oxygen concentrations - see ⑪

In case of controlled ventilation, the inspiratory oxygen concentration depends on frequency and tidal volume. The oxygen concentration which is normally desirable (of 30-50%) can usually be obtained without oxygen reservoir. If higher concentrations are indicated, the oxygen reservoir must be connected to the bag's inlet valve.

Examples of oxygen concentration that can be obtained using an oxygen reservoir at different volumes and frequencies, are shown in ⑪).

A simplified table is printed in the oxygen reservoir bag for fast guidance (see ⑩). The oxygen concentrations indicated are obtained as a minimum as long as the minute ventilation is equal to or less than $300 \times 20 = 6$ l/min.

Example:

Administer 2 l/min of oxygen through the inlet of the oxygen reservoir bag. With a tidal volume of max. 300 ml and a compression frequency of max. 20 per minute, a child receives an oxygen concentration of at least 40%.

6. Description

The Ambu Mark IV Baby resuscitator is based on the original Ambu double wall principle which transfers the ventilation pressure directly to operator's fingers. During normal operation, the outer cover expands, and the airway pressure is thereby limited to approx. 4.5 kPa (45 cmH₂O).

The Ambu Mark IV Baby consists of the following main parts: resuscitator bag 12.1, patient valve 12.2, inlet valve 12.3, reservoir tube 12.4 or closed reservoir 12.5 with adapter. Mask, tubes and PEEP 10 valve are accessories, see ⑫.

6.1. Bag - see ⑬

The Ambu Mark IV Baby resuscitator is designed in accordance with the original Ambu double wall principle. The resuscitator bag consists of a self-expanding inner bag made from shaped silicone rubber with a thin-walled outer cover, which fits airtight around the outlet connectors.

The outer cover is made from highly resilient silicone rubber which expands in operation and thus reduces the risk of creating dangerously high pressures in the patient's airways. The bag is provided with a hanging strap with nameplate (see ⑬). The plate can be stamped or engraved.

13.1: Connector, 13.2: Outer cover, 13.3: Inner bag, 13.4: Inlet valve, 13.5: Hanging strap incl. a window with a snap fastener, 13.6: Nameplate, 13.7: Outlet connector.

6.2. Patient valve - see ⑭

The function of the valve is described in Section 4.

⑭ 14.1: 30 mm expiratory male connector, 14.2: valve disc, 14.3: 22/15 mm patient connector, 14.4: manometer port

6.3. Inlet valve - see ⑮

The inlet valve has an oxygen tube connector and connection for standard 22 mm open reservoir tube.

⑮ 15.1: Oxygen tube connector for tubes with inside diameter of 6 mm, 15.1: To be closed by adapter when using oxygen reservoir bag, 15.2: connection for reservoir tube.

6.4. Masks - see ⑯

The Ambu Baby resuscitator can be supplied with one or more of the masks below, as desired.

Ambu mask 0A with transparent shell (16.1) and removable, extremely soft mask bead made from rubber with a sealed foam inner lining. Connection cone: 15 mm ISO. Dead volume: approx. 25 ml.

Rendell-Baker mask in black rubber (16.2). Connection for 22 mm ISO patient connector.

Size: 0: Dead volume: approx. 4 ml

Size: 1: Dead volume: approx. 6 ml

Other mask types can be supplied as accessories.

6.5. Oxygen reservoir bag (optional) - see ⑰

The oxygen reservoir bag (17.12) is connected to a reservoir housing (17.10) with cap (17.11), which fits the inlet valve connector on the resuscitator using an adapter (17.9).

The valve has two silicone rubber valve discs that enable intake of external air when the bag is empty, and discharge of excess oxygen when the bag is full.

7. Cleaning - disinfection - sterilisation

⑰ shows the resuscitator disassembled for cleaning, disinfection and sterilisation purposes.

Note: The parts should not be further disassembled than shown in the picture.

Clean - disinfect - sterilise after each patient or in accordance with local procedures.

⑰ 17.1: Patient valve, 17.2: OA mask, 17.3: Rendell-Baker mask, 17.4: Inner lining, outer cover and connector, 17.5: Inlet valve, 17.6: Reservoir tube, 17.7: Oxygen reservoir bag and adapter, 17.8: Oxygen supply tube.

Clean - disinfect - sterilise the entire resuscitator equipment after use on patients or in environments containing infectious diseases.

Applicable methods	Cleaning		Disinfection - sterilisation				
	Washing		Disinfection			Autoclaving	
● Applicable ○ Not applicable	Manual washing	Washing machine	Heat disinfection	Boiling	Chemical	121°C (250°F)	134°C (273°F)
Parts							
17.1	●	●	●	●	●	●	●
17.2	●	○	○	○	●	○	○
17.3	●	●	●	●	●	●	●
17.4	●	●	●	●	●	●	●
17.5	●	●	●	●	●	●	●
17.6	○	○	○	○	●	○	○
17.7	●	●	●	●	●	●	●
17.8	●	●	●	●	●	●	●

Detergents and chemical disinfectants

There are a wide range of detergents and chemical disinfectants available on the market. Use only detergents or chemical disinfectants that are compatible with the materials of the resuscitator. Otherwise, the life of the equipment may be reduced. See the materials list in Section 3. Follow the manufacturer's guidelines for dilution and exposure time of the detergent or the chemical disinfectant. Agents containing phenol should be avoided.

On request, information on the compatibility of a selection of detergents and chemical disinfectants with the materials of the resuscitator is provided by the nearest Ambus dealer or Ambu A/S.

Note: Avoid cleaning the product with agents containing phenol. Phenol causes rapid wear and deterioration of the materials or a reduced life of the equipment.

Note: Immediately remove all detergent or disinfectant residues from the resuscitator. These residues may result in rapid wear or a reduced life of the device.

Cleaning, disinfection and sterilisation procedure

The following items should be followed. Use the table to select methods that are appropriate for the individual parts.

7.1. Disassembly of the resuscitator

The resuscitator must always be disassembled before cleaning as shown in ⑰.

Note: Do not attempt to disassemble the outlet connector from the inner bag (permanently assembled).

Note: The parts should not be further disassembled than shown in the picture. See ⑰ for further information.

Note: To remove the outer cover of the bag, first unscrew the inlet valve cap. Then roll off the valve housing with the inlet valve disc, see ⑱.

Snap the connector off the inner lining and expand the outer cover at the back, pulling the inner lining out while the bead in front is pulled over the snap cuff, see ⑱.

Removing the strap ⑲:

Disengage the strap at the back of the transparent nameplate cover. Pull the cap over the strap and remove the nameplate. Pull the strap through the slot.

Note: The parts should not be further disassembled than shown.

7.2. Cleaning of parts

Manual washing or automatic washing machine

Wash the parts in water using an appropriate detergent in accordance with the materials list in Section 3. Rinse all parts thoroughly in clean water to remove all detergent residues.

If the exterior faces of the resuscitator bag are to be cleaned and/or disinfected, check that the detergent and/or disinfectant used is compatible with the materials of the resuscitator. Make sure that all residues are removed completely by rinsing with clean water. Otherwise, the life of the materials can be significantly reduced.

An automatic dishwasher with a programme for washing anaesthetics equipment is applicable.

7.3. Disinfection and/or sterilisation

Select heat treatment or chemical disinfection according to local practice for disinfection and the table of applicable methods.

Washing machine - heat disinfection programme: Automatic washing machines designed for medical equipment are usually equipped with heat disinfection programmes.

Boiling: Use clean water. Heat the parts, and boil them for 10 minutes to disinfect them.

Chemical disinfection (e.g. Cidex, chemicals with sodium hypochlorite etc.): Follow the disinfectant manufacturer's guidelines for dilution and exposure time. Rinse all parts thoroughly in clean water to remove all disinfectant residues.

Note: As it is very difficult to thoroughly clean the oxygen reservoir bag, it is recommended not to immerse it in chemical disinfectants.

Autoclaving: Use standard autoclaving equipment set to 121°C or 134°C (250°F or 273°F), respectively, for the individual parts of the resuscitator according to the table above.

7.4. Drying and cooling

Let the parts dry and cool completely before the resuscitator is reassembled.

7.5. Inspection of parts

After cleaning, disinfection and/or sterilisation, all parts must be checked for damage or wear. Replace if necessary. Certain methods can lead to discolouring of rubber parts without affecting the life of these parts. In case of onset of deterioration of the materials, e.g. crack formation, the parts should be replaced.

In case of autoclaving, the oxygen reservoir bag can look slightly wrinkled. This has no impact on its life or function.

7.6. Assembly and functional testing - see 20

Press the inner bag firmly together while expanding the outer cover opening and slide the inner bag into the outer cover. Snap the outer cover bead over the snap cuff at the front of the inner bag. Allow the inner bag to expand and pull the outer cover into its correct position. Check that the outer cover is not inside out (the lot no. on the outer cover must be on the outside). Also check that the outlet connector is intact and positioned correctly. Snap the connector in place at the front. Press the valve housing with the valve discs into the inner bag and screw the inlet valve disc back on.

Patient valve: Before assembling the patient valve, all parts - especially the valve disc - should be checked. Assemble the patient valve in reverse disassembly order.

When inserting the patient valve disc, check that the butt of the pin goes through the hole in the middle of the patient valve housing, and then the expiration connector is screwed in.

Press the patient valve into the bag while turning it.

Perform test for correct function. See Section 8 of these directions for use.

8. Functional test

8.1. Resuscitator

1) Connect an 0,5 l breathing bag to the patient connector. Compress and release the resuscitator several times.

Squeeze and release the resuscitator repeatedly. Check that the breathing bag is filled.

By continued ventilation, the breathing bag must be visibly filled and emptied. See 21.

8.2. Pressure-limiting valve

WARNING

Perform the test below immediately before every use of the resuscitator. Otherwise, the result may be no or reduced airflow.

- 1) Close the pressure-limiting valve with the locking cap and the patient connector with the thumb, while squeezing the resuscitator to check for tightness and proper valve fitting, see .
- 2) Open the pressure-limiting valve by removing the locking cap, and repeat the procedure. The pressure-limiting valve should now be activated, and it should be possible to hear the expiratory flow from the valve.

8.3. Oxygen reservoir system

Reservoir tube:

Set a flow of approx. 5 l/min through the oxygen supply tube. Check that the oxygen flows out at the end of the reservoir tube. If not, check for a blocked oxygen tube.

Oxygen reservoir bag:

Set a flow of approx. 5 l/min, occlude the patient connector on the patient valve and check that the reservoir bag unfolds. Shut the flow. Compress the oxygen reservoir bag and check that the gas leaks unhindered through the holes in the pressure-limiting valve.

9. Service

The Ambu Mark IV Baby resuscitator requires no scheduled maintenance other than regular cleaning, inspection and testing.

It is recommended to thoroughly check all parts after disassembly for cleaning, disinfection or sterilisation.

If any parts are missing or defective, spare parts are available (see the list in Section 11).

Note: The Ambu Mark IV Baby resuscitator must be checked for correct function in accordance with these instructions for use before being released for use again (see Section 8).

10. Storage

The Ambu Mark IV Baby resuscitator should not be exposed to direct sunlight or elevated temperatures during storage.

Never store the resuscitator in close proximity to products containing oil or grease.

Storage temperature: -40°C to 70°C (40°F to 158°F).

Note: Never squeeze the bag too hard on storage.

In case of long-term storage, the resuscitator should be packed away in a cool place away from direct sunlight.

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