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en | Directions for Use of Suction Catheters

INTENDED USE

The suction catheter is a common medical device used to remove the respiratory secretions and keep the respiratory tract open.

Intended for:

1. Inability to cough and for discharging sputum and when patient suffers from difficult breathing such as in coma, critical patient, after anesthesia patient.
2. First aid to asphyxia caused by drowning or amniotic fluid aspiration.

USAGE METHOD

1. Sterilized and ready to use;
2. Examine the function and tube connection before suction operation;
3. Obey the sterile operation rules; the tube, the dressing, the water should be replaced every day to avoid airway pollution;
4. Insert the tube softly, not hardly, reduce the mucosal contact time on the airway, use discontinuous operation method;
5. Supply oxygen again between 2 suction to avoid hypoxemia. The suction should be stopped immediately if Paroxysmal cough or Arrhythmia take place ;
6. The suction cannot last over more than 10-15 seconds

KONTRAINDICATIONS:

1. It is forbidden for patients allergic to PVC materials to use this product.
2. There are no absolute contraindications, but it is forbidden to suck sputum through nasal cavities for a patient with a basal skull fracture.
3. Relative contraindications: Patients with severe depletion of oxygen, and patients with severe arrhythmia.

MATTERS NEEDING ATTENTION:

1. The effective period of sterilization of this product is 5 years and no use of an expired product is allowed.
2. No use of the product is permitted if the package is damaged.
3. The product is disposable, and is strictly forbidden to reuse the product. After being used, the products shall be discarded to avoid pollution. Disposal of the device in accordance with the applicable regional environment regulations.
4. This product contains DEHP, so it should be used with caution by pregnant women and children.
5. Use only by professional medical personnel.
6. When connecting with extensions or connections, please note the funnel size. We recommend step hose connectors size 7/7 mm (Our REF.: I6 6023)

STORAGE AND TRANSPORT

The suction catheter should be stored in clean condition and kept in good ventilation, keep away from corrosive gas; storage temperature range: 2 - 40 °C.

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servoprax GmbH, Am Marienbusch 9, 46485 Wesel, Germany